

OEC SKILL 14-2**Administration with an Auto-injector: Twinject™**


After ensuring that the medication is prescribed for the patient, remove the green cap from the device.




Remove the red cap from the end of the Twinject™.




Assist the patient by holding the Twinject™ with the gray cap against the outside of the patient's thigh, halfway between the hip and the knee, and inject the dose of epinephrine.

OEC SKILL 14-3**Administration with an Auto-injector: Twinject™ Additional Dose**

To administer a second dose using the Twinject™, start by removing the gray cap.



Remove the yellow safety collar from the needle plunger.



Assist the patient by injecting the second dose into the outside of the patient's other thigh.